



 **Healer**™

# ORTHOPEDICS and DIABETIC FOOTWEAR

Exclusively from  
 **I-RUNNER**





Our friendly and knowledgeable Customer Service Representatives are available to take your orders or answer any questions between the hours of 8:00 a.m. and 5:00 p.m. EST.

### To place an order, contact Customer Service

Phone 317-689-7711  
(toll-free) 855-880-7463  
Fax 317-689-7712  
Email [orders@irunners.com](mailto:orders@irunners.com)

### Shipping

Orders are shipped from our centrally located distribution center in Indiana. Orders will ship within 48 hours from the time they are received. All orders are shipped UPS Ground unless specified otherwise.

### Healer Warranty

Healer guarantees all products to be free from defects in materials and workmanship for a period of six months from date of purchase. Defective products will be replaced or credited. Please have the defective product available for return if requested.

### I-Runner Warranty

There is a three-month manufacturers warranty covering defects on all shoes. Defective shoes will be replaced. All other shoes will be evaluated on a case-by-case basis. This warranty is guaranteed for product defect and not wear and tear from everyday use over a period of time.

### Product Return Policy

Returned products will be accepted up to 60 days from date of invoice, as long as they are in their original condition. After 60 days, all returns will have a 15% restocking fee. No returns will be accepted after 120 days.

Returned products should be shipped freight paid to:

Healer Health  
Attn: Returns Dept.  
7271 Mayflower Park  
Zionsville IN 46077

Credit will be applied to the account once the product has been received, inspected and approved for restocking.



## TABLE OF CONTENTS

<i>Cervical</i> .....	4
<i>Shoulders</i> .....	5
<i>Back</i> .....	7
<i>Wrists</i> .....	12
<i>Knees</i> .....	15
<i>Ankles</i> .....	20
<i>Ankle Walkers</i> .....	23
<i>Plantar Fasciitis</i> .....	26
<i>Cast Boots/Post-Op Shoes</i> .....	27
<i>Diabetic/Therapeutic Shoes</i> .....	30
<i>Slip Resistant Shoes</i> .....	32
<i>Heat-Moldable Insoles</i> .....	33



## Cervical Collar

The Cervical Collar provides comfortable support of the cervical spine in a neutral position. It is made of medium density foam with a cotton stockinette cover.

### EXAMPLES OF USE

- Post-operative rehabilitation
- Mild cervical strains

### PRODUCT FEATURES

- Chin cutaway.
- Universal sizing.

Part #	Description
01-9001	Cervical Collar – 2"
01-9002	Cervical Collar – 3"
01-9003	Cervical Collar – 4"



## Clavicle Strap

The Clavicle Strap is a lightweight brace that holds the shoulder in correct alignment to stabilize the clavicle and to help promote healing.

### EXAMPLES OF USE

- Clavicle fractures
- Postural control


### PRODUCT FEATURES

- Contoured foam padding throughout the straps for patient comfort.
- Figure-8 design makes it easy to adjust the brace for a customized fit.

Part #	Description
01-9004-X	Clavicle Strap

X =	Size	Shoulder Circ.*
1	XS	15"-21"
2	SM	19"-28"
3	MD	24"-34"
4	LG	30"-38"
5	XL	34"-40"



Measure from center of shoulder blades, under the arm and around the top of the shoulder.

\*Measurement taken from center of shoulder blades and around shoulder



## Elastic Shoulder Immobilizer

The Elastic Shoulder Immobilizer provides immobilization by holding the arm close to the body at the arm and the wrist.

### EXAMPLES OF USE

- Shoulder injuries or instabilities
- Shoulder dislocations

### PRODUCT FEATURES

- Elastic waist band.
- Foam wrist pad and nylon humeral cuff with contact closure.
- Fits right or left.

Part #	Description
01-8001-X	Elastic Shoulder Immobilizer – Men's
01-8002-X	Elastic Shoulder Immobilizer – Women's

X =	Size	Men's Rib Measurement	Women's Rib Measurement
2	SM	32"-36"	24"-28"
3	MD	36"-40"	28"-32"
4	LG	40"-44"	32"-36"
5	XL	44"-48"	36"-42"

## Shoulder Sling Series

Healer Health's Shoulder Sling Series provides multiple options for post-operative care—immobilization, abduction and external rotation. The innovative design creates maximum comfort and support for the patient recovering from shoulder surgery or injury.

### EXAMPLES OF USE

- Anterior glenohumeral dislocation/subluxation
- Capsular shifts
- Posterior shoulder repair
- Post-operative Bankart procedures
- Soft tissue strains/repairs
- Rotator cuff repairs

### PRODUCT FEATURES

- Superior fit through contoured forearm design (1).
- Quick-release buckles (2) for easy application and removal—single-handedly!
- Breathability through advanced moisture-wicking materials.
- Premium shoulder strap padding (3) for increased comfort.
- Universal left or right.
- Standard and Neutral pillow options—includes exercise ball (4).

Part #	Description
01-8003-X	Shoulder Abduction Sling, Std
01-8004-X	Shoulder Abduction Sling, Neutral
01-8005-X	Shoulder Sling, Premium
01-8006-X	Shoulder Immobilizer, Premium

X =	Size	Length from Elbow*
2	SM	Up to 11"
3	MD	11"-13"
4	LG	13"-15"
5	XL	15"+

\*Length from elbow to base of index finger



Shoulder Abduction Sling, Std

Shoulder Abduction Sling, Neutral



Shoulder Sling, Premium

Shoulder Immobilizer, Premium



Standard Pillow




Neutral Pillow



1



2



3



4

## Shoulder Immobilizer – Standard or Deluxe

The Shoulder Immobilizer comes in Standard and Deluxe options and provides patient comfort, while immobilizing the arm and shoulder.

### EXAMPLES OF USE

- Shoulder injuries or instabilities

### PRODUCT FEATURES

- Soft cotton/poly blend for increased comfort.
- Removable foam waist strap (standard only).
- D-ring on strap for easy adjustment.
- Fits right or left arm.
- Deluxe option has fully padded shoulder strap.



Part #	Description
01-8007-X	Shoulder Immobilizer – Standard
01-8008-X	Shoulder Immobilizer – Deluxe

X =	Size	Length
1	PED	10.5"
2	SM	12.5"
3	MD	15.0"
4	LG	18.0"
5	XL	19.5"

## Deluxe Sling and Swathe

The Deluxe Sling and Swathe safely positions the arm close to the body for immobilization.

### EXAMPLES OF USE

- Shoulder injuries or instabilities
- Shoulder dislocations

### PRODUCT FEATURES

- Cool, lightweight pile/foam construction provides patient comfort during rehab.
- Two over-the-shoulder crisscross straps attach to the sling with contact closures.
- One size fits up to 54" waist.



Part #	Description
01-8009	Deluxe Sling and Swathe



## Envelope Arm Sling

The Envelope Arm Sling provides support to injured arms following trauma or surgery.

### EXAMPLES OF USE

- Arm injuries or instabilities
- Post-surgery

### PRODUCT FEATURES

- Soft cotton/poly blend for increased comfort.
- Contact closure and D-rings for easy adjustment.
- Fits right or left arm.
- Deluxe option has fully padded shoulder strap.

Part #	Description	
01-8010-X	Envelope Arm Sling – Standard	
X =	Size	Length
1	PED	10.0"
2	SM	13.0"
3	MD	14.5"
4	LG	17.0"
5	XL	19.5"



## Mesh Back Support

The Mesh Back Support provides compression for lumbar stabilization.

### EXAMPLES OF USE

- Low back pain
- Low back sprains and strains
- Low back support and stabilization

### PRODUCT FEATURES

- 8" low-profile brace.
- Constructed of breathable, mesh material.
- Molded lumbar compression pad.
- Double side pulls for added support.

Part #	Description	
01-7001-X	Mesh Back Support	
X =	Size	Waist Circ.
2	SM	26"-30"
3	MD	30"-34"
4	LG	34"-38"
5	XL	38"-42"
6	2XL	42"-50"
7	3XL	50"-56"



## Elastic Back Support – Black

The Elastic Back Support provides compression and excellent support to the lower back.

### EXAMPLES OF USE

- Low back pain
- Low back sprains and strains
- Low back support and stabilization

### PRODUCT FEATURES

- 8" low-profile brace.
- Constructed of durable elastic, with contact closures.
- Double side pulls for added support.

Part #	Description	
01-7002-X	Elastic Back Support	
X =	Size	Waist Circ.
2	SM	26"-30"
3	MD	30"-34"
4	LG	34"-38"
5	XL	38"-42"

## Introducing the Healer BackForce Fitting System

The BackForce Fitting System is a multifaceted modular bracing system that provides sizing options for the majority of patients. The complete BackForce Fitting System provides four universal sized medicare-covered back braces that replace twenty skus in your stockroom.

### Why BackForce?

The modularity of the BackForce provides the ability to go with a variety of components for the level and support needed.

Sizing ranges from 25" to 55" bringing added value with inventory control.

### One Brace

Fits waists from 25" - 55"

No more having to order Small, Medium, Large, XL or 2X. No more having to worry about a Universal that is big and bulky from excess material.

### Not A Typical Universal

What makes the BackForce Fitting System so unique is it's adjustable back strap.

Adjusting the size of the belt is as simple as sliding the lateral extensions to the desired size.

You remove the hook and loop barriers, slide the belt extensions and secure them again in the new position with the fasteners.

## Now you can take back control of your inventory



Victory Brace  
HCPCS Lo626



Freedom Brace  
HCPCS Lo627



Inspire Brace  
HCPCS Lo631



Encourage Brace  
HCPCS Lo637

**Universally Sized & Medicare Billable**



## Healer's Dynamic Pulley Closure

There are some other pulley type braces on the market but the BackFource Fitting System has the lowest pull force threshold, meaning it is much easier to adjust the compression of our brace than some of the others. If you have reduced muscle strength you will certainly notice how much easier these braces are to adjust.

It's so easy to adjust and can be fine-tuned, with just one hand throughout the course of your day- loosen it while driving, tighten it when walking; you never have to take the brace off once you put it on.

The advantage of a back support with our proprietary pulley system is that it allows the individual to adjust their level of support throughout the day within a second with a simple pull.

**One-pull adjustment**



**Low friction pulley**



## Indications:

- Post-operative stabilization protocol following spinal fusion, laminectomy/laminotomy, foraminotomy, laproscopic disk replacement, IDET procedures, multi-level decompression.
- Control of unintended thoracic flexion posture.
- Chronic & mechanical low back pain.
- Thoracic mechanical back pain..
- Promotion of spinal alignment to aid in pain reduction associated with osteoarthritis and other spinal anomalies.
- Early onset through variable stages of osteoporosis.
- Compression fractures.
- Facet syndrome.
- Thoracolumbar injury.
- Degenerative disc



## Victory

Suggested HCPCS Lo626

Type: LO

Style: Posterior panel

Purpose: Stability from L1 to S1

Part #	Description	Waist Size
01-0626-16	Victory Back Brace System	25" - 55"



## Freedom

Suggested HCPCS Lo627

Type: LO

Style: Rigid anterior and posterior panels

Purpose: Stabilization from L1 to S1

Part #	Description	Waist Size
01-0627-16	Freedom Back Brace System	25" - 55"



## Encourage

Suggested HCPCS Lo637

Type: LSO

Style: Rigid anterior and posterior panels

Purpose: Stabilization from T9 to S1 with lateral control

Part #	Description	Waist Size
01-0637-16	Encourage Back Brace System	25" - 55"



## Inspire

Suggested HCPCS Lo631

Type: LSO

Style: Rigid anterior and posterior panels

Purpose: Stabilization from T9 to S1

Part #	Description	Waist Size
01-0631-16	Inspire Back Brace System	25" - 55"





## Tennis Elbow

The Tennis Elbow strap provides targeted compression to relieve pain caused by tennis elbow.

### EXAMPLES OF USE

- Tendonitis
- Lateral/medial epicondylitis

### PRODUCT FEATURES

- 2" wide nylon strap with D-ring closure.
- Universal sizing.

Part #	Description	Size
01-6001	Tennis Elbow	Univ



## WristGUARD 6", 8", 10" Wrist Brace

The WristGUARD Wrist Brace is a great solution for wrist support and immobilization. Constructed of a lightweight, durable material, this brace offers maximum comfort.

### EXAMPLES OF USE

- Wrist immobilization
- Sprains and strains
- Tendonitis
- Post-cast healing

### PRODUCT FEATURES

- Durable, lightweight material is Lycra-lined for breathability and comfort.
- Removable and malleable aluminum palmar stay provides correct fit and proper support.
- Easy-to-use contact closures.
- Available in three lengths—6", 8" and 10".

Part #	Description
01-6002-X	WristGUARD Wrist Brace, Right – 6"
01-6003-X	WristGUARD Wrist Brace, Left – 6"
X =	Size      Wrist Circ.
2	SM      6.25"-7.25"
3	MD      7.25"-8.00"
4	LG      8.00"-9.00"
5	XL      9.00"-10.00"

Part #	Description
01-6004-X	WristGUARD Wrist Brace, Right – 8"
01-6005-X	WristGUARD Wrist Brace, Left – 8"

X =	Size	Wrist Circ.
1	XS	5.25"-6.25"
2	SM	6.25"-7.25"
3	MD	7.25"-8.00"
4	LG	8.00"-9.00"
5	XL	9.00"-10.00"

Part #	Description
01-6006-X	WristGUARD Wrist Brace, Right – 10"
01-6007-X	WristGUARD Wrist Brace, Left – 10"

X =	Size	Wrist Circ.
1	XS	5.25"-6.25"
2	SM	6.25"-7.25"
3	MD	7.25"-8.00"
4	LG	8.00"-9.00"
5	XL	9.00"-10.00"



## Wrist Lacer

The Wrist Lacer is a lightweight, low profile brace that immobilizes the wrist, while allowing full function of the hand. The dual lacing closure system, adjustable dorsal stay and malleable palmar stay provide a customized fit and superior support.

### EXAMPLES OF USE

- Wrist immobilization
- Sprains and strains
- Tendonitis
- Post-cast healing

### PRODUCT FEATURES

- Constructed of a soft, comfortable foam material.
- Adjustable dorsal stay.
- Malleable and removable palmar stay.
- Dual pull lacing closure for easy application.
- Universal sizing.

Part #	Description	Size
01-6018	Wrist Lacer – Right	Univ
01-6019	Wrist Lacer – Left	Univ



## Universal 8" Wrist Brace

The Universal 8" Wrist Brace is constructed of a durable, lightweight material that provides both support and comfort. The universal sizing reduces inventory needs.

### EXAMPLES OF USE

- Wrist immobilization
- Sprains and strains
- Tendonitis
- Post-cast healing

### PRODUCT FEATURES

- Lycra-lined for breathability and patient comfort.
- Removable and malleable aluminum palmar stay provides correct fit and proper support.
- Universal sizing.

Part #	Description
01-6010	Universal 8" Wrist Brace – Right
01-6011	Universal 8" Wrist Brace – Left



## Wrist Lacer with Thumb Spica

The Wrist Lacer with Thumb Spica is a lightweight, low profile brace that offers excellent support and immobilization. The dual lacing closure system, adjustable dorsal stay and malleable palmar stay provide a customized fit and superior support.

### EXAMPLES OF USE

- Wrist and thumb sprains and strains
- Wrist immobilization
- DeQuervain's syndrome
- Post-cast healing

### PRODUCT FEATURES

- Constructed of a soft, comfortable foam material.
- Adjustable dorsal stay.
- Malleable and removable palmar stay for customized fit.
- Dual pull lacing closure for easy application.
- Adjustable thumb strap for patient comfort.
- Universal sizing.

Part #	Description	Size
01-6012	Wrist Lacer with Thumb Spica – Right	Univ
01-6013	Wrist Lacer with Thumb Spica – Left	Univ



## WristGUARD 8" Thumb Spica

The WristGUARD 8" Thumb Spica is a great solution for wrist/thumb support and immobilization. Constructed of a lightweight, durable material, this brace offers maximum comfort.

### EXAMPLES OF USE

- Wrist and thumb sprains and strains
- Wrist immobilization
- DeQuervain's syndrome
- Post-cast healing

### PRODUCT FEATURES

- Durable, lightweight material is Lycra-lined for breathability and comfort.
- Removable and malleable aluminum palmar stay provides support for palmar surface of wrist and extensor surface of thumb.
- Easy-to-use contact closures.

Part #	Description
01-6014-X	WristGUARD 8" Thumb Spica – Right
01-6015-X	WristGUARD 8" Thumb Spica – Left

X =	Size	Wrist Circ.
1	XS	5.25"–6.25"
2	SM	6.25"–7.25"
3	MD	7.25"–8.00"
4	LG	8.00"–9.00"
5	XL	9.00"–10.00"



## Universal 8" Thumb Spica

The Universal 8" Thumb Spica is constructed of a durable, lightweight material that provides both support and comfort. The universal sizing reduces inventory needs.

### EXAMPLES OF USE

- Wrist and thumb sprains and strains
- Wrist immobilization
- DeQuervain's syndrome
- Post-cast healing

### PRODUCT FEATURES

- Lycra-lined for breathability and patient comfort.
- Removable and malleable aluminum palmar stay provides correct fit and proper support.
- Adjustable dorsal stay.
- Universal sizing.

Part #	Description	Size
01-6016	Universal 8" Thumb Spica – Right	Univ
01-6017	Universal 8" Thumb Spica – Left	Univ

## Restore – Anterior Frame



Healer’s Restore is a feature-rich Functional Knee Brace that creates superior fit and comfort, resulting in increased patient compliance. Ideal for all activity levels.

### EXAMPLES OF USE

- Moderate to severe ACL instabilities
- Prophylactic/bilateral use
- MCL/LCL instabilities
- ACL reconstructions

### PRODUCT FEATURES

- Opti-Fit technology (1), located on the medial thigh cuff, creates an intimate fit and reduces migration by allowing the natural motion of the knee during flexion and extension.
- Quick-release buckles (2) are internally mounted for easy on and off.
- Straps pivot to create a contoured fit and have gecko tabs (3) with large embossed numbers, making the brace easy to apply.
- Aircraft grade aluminum provides superior strength, yet is extremely lightweight (weighs an average of only 18.5 oz.).
- ROM hinge with flexion and extension control (stop kit included).
- One year warranty on frames and hinges.



## Restore – Step Thru



Part #	Description
01-5001-X	Restore, Ant Frame, Std, OTS, Right
01-5002-X	Restore, Ant Frame, Std, OTS, Left
01-5003-X	Restore, Step-thru, Std, OTS, Right
01-5004-X	Restore, Step-thru, Std, OTS, Left
01-5005-X	Restore PCL, Ant Frame, Std, OTS, Right
01-5006-X	Restore PCL, Ant Frame, Std, OTS, Left

X =	Size	Circ. Above*	Circ. Below**
1	XS	13.00"-15.50"	10.00"-11.75"
2	SM	15.50"-18.00"	11.75"-13.50"
3	MD	18.00"-20.50"	13.50"-15.50"
4	LG	20.50"-23.00"	15.50"-17.50"
5	XL	23.00"-26.00"	17.50"-19.25"
6	2XL	26.00"-29.00"	19.25"-21.50"

\*Circumference measurement taken 6" above mid-patella.

\*\*Circumference measurement taken 6" below mid-patella.

## OA Restore – Medial or Lateral

Utilizing the most innovative materials and design, our lightweight, low profile OA Restore Knee Braces are designed for patients suffering from medial or lateral compartmental osteoarthritis. They feature a Dynamic Unloading Hinge, which provides a varus or valgus force to increase the load in terminal extension. This is accomplished by lengthening the strut on the affected side, changing the angle of the brace and

### PRODUCT FEATURES

- Dynamic Unloading Hinge (1) unloads the knee—and the pain!
- ROM Hinge (2) with adjustable flexion/extension control—Stop Kit not included.
- Cinching Dial (3) makes incremental adjustments a cinch.
- Quick Release Buckles (4) for easy on and off.
- Load Dispersion Cuffs (5) provide increased comfort.

### EXAMPLES OF USE

- Unicompartamental OA-medial and lateral
- Chondral defect procedures
- Meniscal procedures
- Ligament Instabilities

Part #	Description
01-5026-X	OA Restore, OTS, Medial, Right
01-5027-X	OA Restore, OTS, Medial, Left
01-5028-X	OA Restore, OTS, Lateral, Right
01-5029-X	OA Restore, OTS, Lateral, Left

X =	Size	Circ. Above*	Circ. Below**
1	XS	13.00"-15.50"	10.00"-11.75"
2	SM	15.50"-18.00"	11.75"-13.50"
3	MD	18.00"-20.50"	13.50"-15.50"
4	LG	20.50"-23.00"	15.50"-17.50"
5	XL	23.00"-26.00"	17.50"-19.25"
6	2XL	26.00"-29.00"	19.25"-21.50"

\*Circumference measurement taken 6" above mid-patella.

\*\*Circumference measurement taken 6" below mid-patella.



Medial

Lateral

## Healoprene or Neoprene

## Hinged Knee Braces

### PRODUCT FEATURES

Choice - Two materials (Healoprene or Neoprene), two options (Short 13" or Tall 16"). Choose the option that best meets your patient's needs.

Coolness - Our new Healoprene material gives these braces increased breathability and moisture-wicking properties, resulting in a cooler feel.

Comfort - The hinged knee braces utilize a new material in the popliteal area with advanced stretch characteristics. Minimizes bunching, maximizes comfort!

Part #	Description
01-5037-X	Hinged Knee Brace, Short, Healoprene, Wrap
01-5038-X	Hinged Knee Brace, Tall, Healoprene, Wrap
01-5039-X	Hinged Knee Brace, Short, Neo
01-5040-X	Hinged Knee Brace, Tall, Neo, Wrap

XX =	Size	Circ. Above*
1	XS	13.00"-15.50"
2	SM	15.50"-18.00"
3	MD	18.00"-20.50"
4	LG	20.50"-23.00"
5	XL	23.00"-26.00"
6	2XL	26.00"-29.00"

\*Circumference measurement taken 6" above mid-patella.



Hinged Knee Brace,  
Short, Wrap

Hinged Knee Brace,  
Tall, Wrap





Cool Pads

## Rehab Post-Op Knee

The Rehab Post-Op Knee Brace offers a simple-to-use design, which provides protected and controlled range-of-motion for patients recovering from knee surgery.

### EXAMPLES OF USE

- Post-op/post injury immobilization with ROM control

### PRODUCT FEATURES

- Easy-to-use hinge settings (1).
- Drop lock button (2): 0, 15, 30, 45 degrees.
- Adjustable struts (3) from 17" to 28" length.
- Quick-release buckles (4).
- Available in three soft goods styles: Cool Wrap, Cool Pads and Premium Full Foam.

Part #	Description	Circ. Above*
01-5041	Rehab – Cool Wrap, Regular	Up to 27"
01-5042	Rehab – Cool Wrap, XL	27"-35"
01-5043	Rehab – Cool Pads, Regular	Up to 27"
01-5044	Rehab – Cool Pads-XL	27"-35"
01-5045	Rehab – Premium Full Foam, Regular	Up to 27"
01-5046	Rehab – Premium Full Foam, XL	27"-35"

\*Circumference measurement taken 6" above mid-patella.

Cool Wrap

Premium Full Foam



## Knee Sleeve – Economy

### EXAMPLES OF USE

- Mild knee sprains
- Chondromalacia
- Mild patellofemoral osteoarthritis
- Osgood-Schlatter's disease

### PRODUCT FEATURES

- Constructed of neoprene for compression and support.
- Open or closed patella.

Part #	Description
01-5047-X	Knee Sleeve, Open Pat
01-5048-X	Knee Sleeve, Closed Pat

X =	Size	Circ. Above*
2	SM	15.00"-17.50"
3	MD	17.50"-19.50"
4	LG	19.50"-22.50"
5	XL	22.50"-24.50"
6	2XL	24.50"-27.00"

\*Circumference measurement taken 6" above mid-patella.



## Hinged Knee Brace – Economy

### PRODUCT FEATURES

- 1/8" neoprene construction.
- Open patella and closed popliteal.
- Hinges are molded into the brace to maintain a low profile.

Part #	Description	
01-5049-X	Hinged Knee Brace – Economy	
X =	Size	Circ. Above*
2	SM	15.00"-17.50"
3	MD	17.50"-19.50"
4	LG	19.50"-22.50"
5	XL	22.50"-24.50"
6	2XL	24.50"-27.00"

\*Circumference measurement taken 6" above mid-patella.



## Wraparound Hinged Knee Brace – Economy

### PRODUCT FEATURES

- 1/8" neoprene construction.
- Open patella and open popliteal.
- Open popliteal eliminates bunching and increases comfort.
- Front contact closures allow for hassle-free application.
- Wraparound style accommodates all leg shape variations.
- Hinges are removable if desired.

Part #	Description	
01-5050-X	Wraparound Hinged Knee Brace	
X =	Size	Circ. Above*
2	SM	15.50"-18.00"
3	MD	18.00"-20.50"
4	LG	20.50"-23.00"
5	XL	23.00"-26.00"
6	2XL	26.00"-29.00"
7	3XL	29.00"-32.00"
8	4XL	32.00"-35.00"

\*Circumference measurement taken 6" above mid-patella



## Patella Stabilizer – Economy

### PRODUCT FEATURES

- 1/8" neoprene construction.
- Open patella and closed popliteal.
- Molded hardware to maintain a low profile.
- Soft cushioned liner to offer increased comfort.

Part #	Description	
01-5051-X	Patella Stabilizer Universal Buttress	
01-5052-X	Patella Stabilizer Donut Buttress	
X =	Size	Circ. Above*
2	SM	15.00"-17.50"
3	MD	17.50"-19.50"
4	LG	19.50"-22.50"
5	XL	22.50"-24.50"
6	2XL	24.50"-27.00"

\*Circumference measurement taken 6" above mid-patella.



## 3-Panel Knee Immobilizer

The 3-Panel Knee Immobilizer has a three-piece design for a precise fit and is fully adjustable to accommodate varying leg circumferences.

### EXAMPLES OF USE

- Post-operative knee immobilization
- Collateral ligament sprains
- Patella tendon repairs

### PRODUCT FEATURES

- Fully adjustable side panels.
- Lightweight rigid posterior stays, with adjustable medial/lateral stays.
- Universal circumference fits up to 29".

Part #	Description
01-5053	3-Panel Knee Immobilizer – 12"
01-5054	3-Panel Knee Immobilizer – 16"
01-5055	3-Panel Knee Immobilizer – 18"
01-5056	3-Panel Knee Immobilizer – 20"
01-5057	3-Panel Knee Immobilizer – 24"



## Foam Knee Immobilizer

The Foam Knee Immobilizer is designed for quick application and provides universal sizing.

### EXAMPLES OF USE

- Post-operative knee immobilization
- Collateral ligament sprains
- Patella tendon repairs

### PRODUCT FEATURES

- Medial/lateral aluminum stays are removable and moldable for placement and support.
- Posterior aluminum stays are removable and contoured for added rigidity.
- 3" wide compression straps are easy to fasten and provide added support.
- Available in standard and cotton-lined styles.

Part #	Description	Length	Circ. Above*
01-5058	Foam Knee Immobilizer	14"	26"
01-5059	Foam Knee Immobilizer	16"	26"
01-5060	Foam Knee Immobilizer	20"	26"
01-5061	Foam Knee Immobilizer	22"	26"
01-5062	Foam Knee Immobilizer	24"	26"

\*Circumference measurement taken 6" above mid-patella.



## Hinged Ankle Brace

The Hinged Ankle Brace offers support and stability to the ankle and is ideal for transitioning out of a Walking Boot.

### EXAMPLES OF USE

- Sub-acute ankle sprains
- Chronic ankle instabilities
- Ankle injury rehabilitation
- Prophylactic use

### PRODUCT FEATURES

- Hinged side plates and foot plates control inversion/eversion.
- Contoured foot plates provide arch support and stability.
- High density foam pads form to ankle to offer increased comfort.
- Heat moldable shell.
- Fits easily into shoe.

Part #	Description		
01-4003-X	Hinged Ankle Brace – Right		
01-4004-X	Hinged Ankle Brace – Left		
X =	Size	Men's Shoe Size	Women's Shoe Size
1	XS	–	4 – 6
2	SM	4 – 7	7 – 9
3	MD	8 – 10	10 – 12
4	LG	11 – 13	13 – 15



## Foam Stirrup Ankle

The Foam Stirrup Ankle provides stability and support to unstable or injured ankles.


### EXAMPLES OF USE

- Sub-acute ankle sprains
- Chronic ankle instabilities
- Prophylactic use

### PRODUCT FEATURES

- Polypropylene outer shells are anatomically contoured to fit right or left.
- Unique memory foam pads provide uniform compression, a customized fit and maximum comfort.
- A soft, adjustable foot piece connects the shells.

Part #	Description	Height
01-4002-2	Foam Stirrup Ankle – SM	8.75"
01-4002-3	Foam Stirrup Ankle – Trainer	9.50"
01-4002-4	Foam Stirrup Ankle – Regular	10.00"



WHERE TO MEASURE



## Lock-Down Ankle Brace

The Lock-Down Ankle Brace features non-stretch nylon figure-8 straps, which help control abnormal inversion/eversion of the ankle. This low profile brace is ideal for all forms of activity.

### EXAMPLES OF USE

- Sub-acute ankle sprains
- Chronic ankle instabilities
- Prophylactic use

### PRODUCT FEATURES

- Low profile design fits into any shoe.
- Figure-8 nylon straps provide inversion/eversion control.
- Elastic cuff wraps around the ankle for added compression and support.
- Removable medial/lateral stays provide increased stability.
- Fits either right or left ankle.

Part #	Description
01-4005-X	Lock-Down Brace

X =	Size	Ankle Circ.*
1	XS	10.00"-11.00"
2	SM	11.00"-12.00"
3	MD	12.00"-13.00"
4	LG	13.00"-14.00"
5	XL	14.00"-15.00"
6	2XL	15.00"-16.00"

\*Measurement taken diagonally over the base of the heel and across the ankle bone. See diagram.



## Stabilizing Ankle Brace

The Stabilizing Ankle Brace is constructed of non-stretch ballistic nylon for comfort, performance and durability.

### EXAMPLES OF USE

- Sub-acute ankle sprains
- Chronic ankle instabilities
- Prophylactic use

### PRODUCT FEATURES

- Low profile design fits into any athletic or street shoe.
- Figure-8 nylon straps provide inversion/eversion control.
- Removable medial/lateral stays provide increased stability.
- Fits either right or left ankle.

Part #	Description
01-4006-X	Stabilizing Ankle Brace

X =	Size	Ankle Circ.*
1	XS	10.00"-11.00"
2	SM	11.00"-12.00"
3	MD	12.00"-13.00"
4	LG	13.00"-14.00"
5	XL	14.00"-15.00"

\*Measurement taken diagonally over the base of the heel and across the ankle bone. See diagram.



## Lace-Up Ankle Brace

The Lace-Up Ankle Brace is constructed of sturdy canvas for maximum support.

### EXAMPLES OF USE

- Sub-acute ankle sprains
- Chronic ankle instabilities
- Prophylactic use

### PRODUCT FEATURES

- Sturdy canvas design with removable plastic medial/lateral stays.
- Fits into any athletic or street shoe.
- Fits either right or left ankle.

Part #	Description
01-4007-X	Lace-Up Ankle Brace

X =	Size	Ankle Circ.*
1	XS	10.00"-11.00"
2	SM	11.00"-12.00"
3	MD	12.00"-13.00"
4	LG	13.00"-14.00"
5	XL	14.00"-15.00"

\*Measurement taken diagonally over the base of the heel and across the ankle bone. See diagram.



## Air/Gel Stirrup Ankle

The Air/Gel Stirrup Ankle provides stability and support to unstable or injured ankles.

### EXAMPLES OF USE

- Sub-acute ankle sprains
- Chronic ankle instabilities
- Prophylactic use

### PRODUCT FEATURES

- Polypropylene outer shells are anatomically contoured to fit right or left.
- Air/gel bladder, lined with soft foam, provides medial/lateral support and may be used for cold therapy treatment.
- A soft, adjustable foot piece connects the shells.

Part #	Description	Height
01-4008-1	Air/Gel Stirrup Ankle – Ped	6.00"
01-4008-2	Air/Gel Stirrup Ankle – SM	8.75"
01-4008-3	Air/Gel Stirrup Ankle – Trainer	9.50"
01-4008-4	Air/Gel Stirrup Ankle – Regular	10.00"



## Air Stirrup Ankle

The Air Stirrup Ankle provides stability and support to unstable or injured ankles.

### EXAMPLES OF USE

- Sub-acute ankle sprains
- Chronic ankle instabilities
- Prophylactic use

### PRODUCT FEATURES

- Polypropylene outer shells are anatomically contoured to fit right or left.
- Air chambers, lined with soft foam, may be inflated to provide medial/lateral support.
- A soft, adjustable foot piece connects the shells.

Part #	Description	Height
01-4001-2	Air Stirrup Ankle – SM	8.75"
01-4001-3	Air Stirrup Ankle – Trainer	9.50"
01-4001-4	Air Stirrup Ankle – Regular	10.00"

## Healer Air Walker – Standard and Ankle



The Healer Air Walker provides maximum stability and support following injury or post-operative procedures.

### EXAMPLES OF USE

- Acute ankle sprains
- Soft tissue injuries
- Stable lower leg, foot and ankle fractures
- Post-operative use

### PRODUCT FEATURES

- Plastic molded uprights, with steel reinforcement, for increased durability.
- Medial/lateral air bladders offer customized compression to increase comfort and reduce edema.
- Wide footbed provides increased stability, comfort and accommodation for edema and bandages.
- Shock-absorbing insole reduces impact of heel strike when walking.
- Rocker sole helps promote a natural gait to allow continuation of daily activities.

Part #	Description
01-1010-X	Healer Air Walker – Standard
01-1011-X	Healer Air Walker – Ankle

X =	Size	Men's Shoe Size	Women's Shoe Size
1	XS	–	3 – 5
2	SM	5 – 6	5.5 – 7
3	MD	6.5 – 10	7.5 – 11
4	LG	10.5 – 12	11.5+
5	XL	12.5+	–

## Healer Walker – Standard and Ankle



The Healer Walker provides maximum stability and support following injury or post-operative procedures.

### EXAMPLES OF USE

- Acute ankle sprains
- Soft tissue injuries
- Stable lower leg, foot and ankle fractures
- Post-operative use

### PRODUCT FEATURES

- Plastic molded uprights, with steel reinforcement, for increased durability.
- Wide footbed provides increased stability, comfort and accommodation for edema and bandages.
- Shock-absorbing insole reduces impact of heel strike when walking.
- Rocker sole helps promote a natural gait to allow continuation of daily activities.

Part #	Description
01-1014-X	Healer Walker – Standard
01-1015-X	Healer Walker – Ankle

X =	Size	Men's Shoe Size	Women's Shoe Size
1	XS	–	3 – 5
2	SM	5 – 6	5.5 – 7
3	MD	6.5 – 10	7.5 – 11
4	LG	10.5 – 12	11.5+
5	XL	12.5+	–

## Air Walker – Standard and Ankle



The Air Walker provides maximum stability and support following injury or post-operative procedures.

### EXAMPLES OF USE

- Acute ankle sprains
- Soft tissue injuries
- Stable lower leg, foot and ankle fractures
- Post-operative use

### PRODUCT FEATURES

- Fixed, aluminum uprights allow for the frame to be adjusted to conform better to the contour of the leg.
- Medial/lateral air bladders offer customized compression to increase comfort and reduce edema.
- Shock-absorbing insole reduces impact of heel strike when walking.
- Rocker sole helps promote a natural gait to allow continuation of daily activities.

Part #	Description
01-1003-X	Air Walker – Standard
01-1004-X	Air Walker – Ankle

X =	Size	Men's Shoe Size	Women's Shoe Size
1	XS	–	3 – 5
2	SM	5 – 6	5.5 – 7
3	MD	6.5 – 10	7.5 – 11
4	LG	10.5 – 12	11.5+
5	XL	12.5+	–

## Walker Boot – Standard and Ankle

The Walker Boot provides maximum stability and support following injury or post-operative procedures.

### EXAMPLES OF USE

- Acute ankle sprains
- Soft tissue injuries
- Stable lower leg, foot and ankle fractures
- Post-operative use

### PRODUCT FEATURES

- Fixed, aluminum uprights allow for the frame to be adjusted to conform better to the contour of the leg.
- Shock-absorbing insole reduces impact of heel strike when walking.
- Rocker sole helps promote a natural gait to allow continuation of daily activities.

Part #	Description
01-1001-X	Walker Boot – Standard
01-1002-X	Walker Boot – Ankle



X =	Size	Men's Shoe Size	Women's Shoe Size
1	XS	–	3 – 5
2	SM	5 – 6	5.5 – 7
3	MD	6.5 – 10	7.5 – 11
4	LG	10.5 – 12	11.5+
5	XL	12.5+	–





## Deluxe Air Walker

The Deluxe Air Walker has a circumferential, semi-rigid shell that offers maximum support to the leg/ankle following injury or post-operative procedures.

### EXAMPLES OF USE

- Acute ankle sprains
- Soft tissue injuries
- Stable foot and ankle fractures
- Post-operative use

### PRODUCT FEATURES

- Streamlined, semi-rigid shell is lightweight, yet provides maximum support for the leg and ankle.
- Medial/lateral air bladders offer customized compression to increase comfort and reduce edema.
- Rocker sole helps promote a natural gait to allow continuation of daily activities.
- Shock-absorbing insole reduces impact of the heel strike while walking.
- Toe guard for added protection.

Part #	Description
01-1012-X	The Deluxe Air Walker – Standard
01-1013-X	The Deluxe Air Walker – Ankle

X =	Size	Men's Shoe Size	Women's Shoe Size
1	XS	–	4 – 5
2	SM	5 – 6	5.5 – 7
3	MD	6.5 – 10	7.5 – 11
4	LG	10.5 – 12	11.5 – 13
5	XL	12+	13+



## Rom Walker – Standard and Ankle

The ROM Walker provides protected range-of-motion following injury or post-operative procedures.

### EXAMPLES OF USE

- Acute ankle sprains
- Soft tissue injuries
- Stable lower leg, foot and ankle fractures
- Post-operative use following Achilles tendon repair

### PRODUCT FEATURES

- Aluminum uprights allow for the frame to be adjusted to conform better to the contour of the leg.
- ROM settings between 45 degrees plantar-flexion and 30 degrees dorsi-flexion in 7.5 degree increments.
- Rocker sole helps promote a natural gait to allow continuation of daily activities.

Part #	Description
01-1005-X	ROM Walker – Standard
01-1006-X	ROM Walker – Ankle
01-1007-X	Air ROM Walker – Standard
01-1008-X	Air ROM Walker – Ankle

X =	Size	Men's Shoe Size	Women's Shoe Size
2	SM	5 – 6	5.5 – 7
3	MD	6.5 – 10	7.5 – 11
4	LG	10.5 – 12	11.5+
5	XL	12.5+	–



## Pediatric Walker

The Pediatric Walker is designed to treat stable fractures and ankle sprains.

### PRODUCT FEATURES

- Metal uprights.
- Shock-absorbing insole reduces impact of heel strike when walking.
- Rocker sole helps promote a natural gait to allow continuation of daily activities.

Part #	Description	Years of Age
01-1016S	Pediatric Walker – SM/MD	1 – 2 years
01-1016L	Pediatric Walker – LG	2 1/2 – 6 years
01-1016X	Pediatric Walker – XL	6 1/2 – 9 1/2 years



## Diabetic Walker

The Diabetic Walker is designed to provide comfort, stability and protection.

### PRODUCT FEATURES

- Diabetic insole with removable pegs effectively offloads plantar pressure on the foot.
- Rocker sole helps promote a natural gait to allow continuation of daily activities.

Part #	Description	Men's Shoe Size	Women's Shoe Size
01-10091	Diabetic Walker – XS	–	3 – 5
01-10092	Diabetic Walker – SM	5 – 6	5.5 – 7
01-10093	Diabetic Walker – MD	6.5 – 10	7.5 – 11
01-10094	Diabetic Walker – LG	10.5 – 12	11.5+
01-10095	Diabetic Walker – XL	12.5+	–



## PF Wrap

The PF Wrap is designed to minimize the discomfort of plantar fasciitis.

### EXAMPLES OF USE

- Plantar fasciitis

### PRODUCT FEATURES

- Low profile design fits into most shoes.
- Can be worn day and night.

Part #	Description	X =	Size	Ankle Circ.*
01-2004-X	PF Wrap	2	SM	11"
		3	MD	12"
		4	LG	13"

\*Measurement taken around the widest part of the ankle



## Silicone Heel Cups

### PRODUCT FEATURES

- Made of molded medical-grade soft silicone.
- Helps reduce shock and vibration loads to joints.
- Soft center in heel region helps alleviate pain in sensitive area.

Part #	Description	Size
01-20031	Silicone Heel Cups	Child/Youth
01-20032	Silicone Heel Cups	Adult REG
01-20034	Silicone Heel Cups	Adult LG



## Dorsal Night Splint

The Dorsal Night Splint is a low profile solution for plantar fasciitis treatment.

### EXAMPLES OF USE

- Plantar fasciitis
- Achilles tendonitis
- Drop foot

### PRODUCT FEATURES

- Lightweight, low profile design for increased comfort and compliance.
- Slips on and adjusts easily.
- Non-slip pad on plantar pad for limited ambulation.

Part #	Description
01-2002-X	Dorsal Night Splint

X =	Size	Men's Shoe Size	Women's Shoe Size
2	SM/MD	5 – 9	6 – 10
4	LG/XL	9.5 – 14	10.5 – 13



## Adjustable Night Splint

The Adjustable Night Splint features adjustable straps for controlled application of a prolonged stretch. Non-skid tread for limited ambulation.

### EXAMPLES OF USE

- Plantar fasciitis
- Achilles tendonitis

### PRODUCT FEATURES

- Adjustable flexion straps provide a gentle stretching of the plantar fascia and Achilles tendon.
- Adjusts from 10° to 90° dorsi-flexion.
- Lightweight, sturdy, one-piece flexible shell.
- Cool foam liner, padded calf and ankle straps for additional comfort.

Part #	Description
01-2001-X	Adjustable Night Splint

X =	Size	Men's Shoe Size	Women's Shoe Size
2	Small	Up to 7	Up to 7.5
3	Medium	7.5 – 9	8 – 9.5
4	Large	9.5 – 11.5	10 – 12.5



## Square Toe Post-Op Shoe

The Square Toe Post-Op Shoe provides increased protection and comfort following post-operative procedures and forefoot injuries.

### PRODUCT FEATURES

- Square toe design provides better left/right fit.
- Forefoot closure offers easy on/off and accommodates large dressings.
- Semi-rigid sole.
- Rocker sole.

Part #	Description
01-3002-X	Square Toe Post-Op Shoe

X =	Size	Shoe Size
1	Ped	12 – 1
2	Men's SM	6 – 8
3	Men's MD	8.5 – 10
4	Men's LG	10.5 – 12
5	Men's XL	12.5 – 14
12	Women's SM	4 – 6
13	Women's MD	6.5 – 8
14	Women's LG	8.5 – 10

## Mesh Top Post-Op Shoe



### PRODUCT FEATURES

- Breathable mesh fabric upper with padded heel collar.
- D-ring with loop/lock closure and extended straps accommodate various dressing widths.
- Semi-rigid sole.

Part #	Description
01-3001-X	Mesh Top Post-Op Shoe

X =	Size	Shoe Size
1	Ped	12 – 1
2	Men's SM	6 – 8
3	Men's MD	8.5 – 10
4	Men's LG	10.5 – 12
5	Men's XL	12.5 – 14
6	Men's XXL	14.5+
7	Women's SM	4 – 6
8	Women's MD	6.5 – 8
9	Women's LG	8.5 – 10



## Rigid Sole Post-Op Shoe

### PRODUCT FEATURES

- Breathable mesh fabric upper.
- D-ring with loop/lock closure and extended straps accommodate various dressing widths.
- Rigid sole.

Part #	Description
01-3003-X	Rigid Sole Post-Op Shoe

X =	Size	Shoe Size
1	Ped	12 – 1
2	Men's SM	6 – 8
3	Men's MD	8.5 – 10
4	Men's LG	10.5 – 12
5	Men's XL	12.5 – 14
12	Women's SM	4 – 6
13	Women's MD	6.5 – 8
14	Women's LG	8.5 – 10



## Blue Nylon Post-Op Shoe

### PRODUCT FEATURES

- Lightweight blue nylon upper.
- Extended straps accommodate various dressing widths.
- Rigid sole.

Part #	Description
01-3004-X	Blue Nylon Post-Op Shoe

X =	Size	Shoe Size
1	Ped	12 – 1
2	Men's SM	6 – 8
3	Men's MD	8.5 – 10
4	Men's LG	10.5 – 12
5	Men's XL	12.5 – 14
12	Women's SM	4 – 6
13	Women's MD	6.5 – 8
14	Women's LG	8.5 – 10



## Canvas Cast Shoe

### PRODUCT FEATURES

- Sandal style with canvas upper.
- Open toe and open heel.
- Forefoot and heel loop/lock closure.
- E.V.A. rocker sole.

Part #	Description
01-3005-X	Canvas Cast Shoe

X =	Size	Sole Length
00	Infant	5.25"
01	PED	6.50"
1	XS	7.75"
2	SM	9.50"
3	MD	10.50"
4	LG	11.75"
5	XL	13.50"



## Canvas Cast Boot

### PRODUCT FEATURES

- Sandal style with canvas upper.
- Closed toe and closed heel.
- Loop/lock contact closure.
- E.V.A. rocker sole.

Part #	Description
01-3006-X	Canvas Cast Boot

X =	Size	Sole Length
01	PED	8.00"
1	XS	9.00"
2	SM	10.75"
3	MD	11.50"
4	LG	12.75"
5	XL	13.75"



## Rocker Cast Shoe

### PRODUCT FEATURES

- Canvas upper.
- Open toe and open heel.
- Adjustable forefoot and heel loop/lock closures.
- Rocker sole.

Part #	Description
01-3007-X	Rocker Cast Shoe

X =	Size	Sole Length
1	XS	7.25"
2	SM	9.50"
3	MD	11.00"
4	LG	12.75"
5	XL	13.75"



## Equinus Cast Shoe

### PRODUCT FEATURES

- Sandal style with canvas upper.
- Open toe and open heel.
- Forefoot and heel loop/lock closure.
- E.V.A. rocker sole.

Part #	Description
01-3008-X	Equinus Cast Shoe

XX =	Size	Sole Length
01	PED	6.50"
1	XS	8.00"
2	SM	10.00"
3	MD	10.50"
4	LG	11.50"
5	XL	13.50"



## Vinyl Cast Boot

### PRODUCT FEATURES

- Excellent for use with fiberglass casts or bulky compression bandages.
- Open toe and closed heel.
- Weatherproof, stretchable upper.
- Loop/lock contact closure.
- Rocker sole is designed to help promote a natural gait.

Part #	Description
01-3009-X	Vinyl Cast Boot

XX =	Size	Men's Shoe Size	Women's Shoe Size
1	XS		3 - 5
2	SM		5.5 - 8
3	MD	6 - 8	8.5 - 10
4	LG	8.5 - 11	
5	XL	11.5 - 14	

## I-Runner Men's Comfort Shoes

Treat your feet to comfort with these athletic shoes by I-Runner. Ideal for casual walking, working out or just lounging around the house, every step you take will be ultra-comfortable thanks to the extra-depth design and wide toe box.



## PDAC approved A5500

- Lightweight
- Breathable mesh upper
- Wider seam-free toe box to promote circulation
- Shock absorbing midsole
- Composite shank for added stability
- Additional padding for luxurious comfort and secure fit
- Extra depth and removable foot bed to accommodate inserts
- Molded ankle counter for prevention of pronation and supination
- Slight rocker sole to encourage natural heel to toe gait

Available in sizes 7.5 - 12, 13, 14, 15

## I-Runner Women's Comfort Shoes

The perfect, wear-anywhere shoe, the I-Runner Comfort Line has a flattering style that looks great on everyone. These lightweight and comfortable athletic shoes have a look that transition easily from sports to casual comfort.



**Sophia**  
Black/Pink

**Eliza**  
Grey/Blue

**Miya**  
White/Pink

**Maria**  
Dark Grey/Salmon

### PDAC approved A5500

- Lightweight
- Breathable mesh upper
- Wider seam-free toe box to promote circulation
- Shock absorbing midsole
- Composite shank for added stability
- Additional padding for luxurious comfort and secure fit
- Extra depth and removable foot bed to accommodate inserts
- Molded ankle counter for prevention of pronation and supination
- Slight rocker sole to encourage natural heel to toe gait

Available in sizes 6 - 11

## I-Runner Non-Slip Pro Series

Crafted from full-grain leather or a leather and mesh combo, the Pro Series from I-Runner will allow you to not only walk in comfort but safety as well. The state-of-the-art anti-slip sole provides traction on a variety of surface types. Tested and qualifies as work wear but comfortable enough to wear them everyday. Protect against slips, trips and falls with these non-slip diabetic shoes.

- All of the benefits of the I-Runner Comfort Series plus an extraordinary level of oil and skid resistance for safety in the workplace and at home.
- Perfect for employees in healthcare, food service, supermarketing, factories, fire departments, law enforcement and other service industries.
- Extra depth and removable foot bed to accommodate inserts
- Wider seam-free toe box to promote circulation
- Composite shank for added stability



**Mesh/Leather**  
*Pro Series Slip Resistant*

**Leather**  
*Pro Series Slip Resistant*

## PDAC approved A5500

Men's sizing  
Available in sizes 7.5 - 12, 13, 14, 15

Women's sizing  
Available in sizes 6 - 11





Medicare Covered  
**Cloud Insole**  
PDAC A5512

- The base layer is a proprietary heat-moldable polyurethane, developed specifically for this product in collaboration with a US based chemical company.
- The top layer is a closed cell crosslinked polyethylene foam, also manufactured in the United States.
- Ideal for diabetics.
- Perfect for anyone that wants the ultimate comfort insole.
- Micro Cell top cover makes it a great insole for walking.
- Conforms to the shape of your feet for amazing pressure relief.

Men's sizing

Available in sizes 7.5 - 12, 13, 14, 15

Women's sizing

Available in sizes 6 - 11



Cash 'n Carry  
**Comfort Insole**

- Shock, impact and vibration absorbing inlay
- Top cover prevents bacteria, mold and is non-abrasive
- Durable EVA middle layer provides superior cushioning
- Medial Arch fill supports over pronated flat feet and supinated high arches
- EVA bottom layer shell provides support and cushioning
- High load bearing insert
- Perfect for every one to wear everyday.

Men's sizing

Available in sizes 7.5 - 12, 13, 14, 15

Women's sizing

Available in sizes 6 - 11



***Supporting  
Your Life...***

Physically, mentally, emotionally





1-855-880-7463 | [HEALERHEALTH.COM](http://HEALERHEALTH.COM)

HEALER HEALTH

7271 MAYFLOWER PARK | ZIONSVILLE, IN 46077 | P 317-689-7711 | F 317-689-7712

HEALER PRODUCTS ARE DISTRIBUTED BY I-RUNNER